



# Veterinary Flexible Endoscopy System

# Veterinary Image Processor



## Full-HD Resolution

- Portable touch screen monitor with full-HD (1920×1080p) resolution offering excellent imaging on 13.3" screen



## Plug and Play

- See a live image within seconds after turning it on.
- Long unplugged battery time supports free move.



## Portable

- Easy to carry around through it's ergonomic handle.
- Easy mounted on Scivita Medical's trolley without screw or bracket.
- Easy data transfer with USB, support Type A and Type C.



## Hygienic System

- Easy to clean



## Multi-Function

- External Storage
- Preference Setting
- Operation Setting
- Live View
- File Explorer



## Compatibility

- Compatible with all Scivita Medical's single use videoscope product family.
- Compatible with displays of main stream brands in high quality to display live image on external screen.
- Photo and video recording.

Veterinary Image Processor			
Mould	SVET-300A		
Size(mm)	332×240×50	IPX	Front Side: IP66, Back Side: IP30
Screen Size(mm)	13.3" TFT-LCD	Weight (kg)	2.3
Resolution	1920×1080p	USB Interface	Type-A×1, Type-C×1
Output	HDMI×1, SDI×1	Storage	Max 128GB* + External USB Storage
Special Features	Sharpness Denoise Gain	Battery Operation	4H



## SCIVITA MEDICAL TECHNOLOGY CO.,LTD.

+86-512-8187 7788

www.scivitamedical.com

info@scivitamedical.com

No.2, Qingqiu Street, Suzhou

Industrial Park, Suzhou, Jiangsu

Province, P.R.China, 215127

